

Mimitakara[®]



@goodHearing A1

User Guide


1. ACCESSORIES	P 1
2. FEATURES	P 2
3. SET UP	P 3
3.1) Wearing Instructions	P 3
3.2) Recharging Instruction	P 5
3.3) Switching Between Modes	P 7
4. PAIRING BLUETOOTH	P 9
4.1) Download App.....	P 9
4.2) Bluetooth Connect	P 9
4.3) Notes	P 12
5. SAFETY & MAINTENANCE	P 14
5.1) Earbud Maintenance	P 14
5.2) Battery Maintenance	P 15
5.3) LED Statuses	P 16
6. FAQ	P 18
7. SPECIFICATION	P 21
8. WARNING AND REMINDER	P 22
9. LEGAL INFORMATION	P 23

1. ACCESSORIES

- ① Main device with replaceable earbuds (M size)
- ② 3 x Additional pair of L, S & Mini size earbuds
- ③ Cleaning Brush
- ④ USB Cable



2. FEATURES



- ① Cable Clips
- ② Earphones
- ③ Charging Port
- ④ LED Indicator
- ⑤ Power Switch
- ⑥ Dual Microphone
- ⑦ Magnetic Earphone Holders
- ⑧ Left Volume Controller
- ⑨ Right Volume Controller
- ⑩ Mode / Bluetooth Connect

3. SET UP

3.1) Wearing Instructions

- 1) Simply place the device around your neck, under your shirt collar or on the neck and place the earbuds into your ear.



- 2) The power switch is on the left side of the device facing your neck.



3) The **goodHearing A1** has adjustable volume settings. To adjust the volume of the device, the right and left sides **+** / **-** buttons situated on the outer sides base of the device, increase and decrease the volume respectively to suit your preferred specification

CAUTION: LISTENING AT A HIGH VOLUMES MAY DAMAGE YOUR HEARING

4) We suggest placing the device around your neck, under your collar to avoid it from rubbing against objects in order to reduce static sounds (Ensure the microphone is not obstructed for optimum sound quality)

3.2) Recharging Instruction

- 1) Charging the battery before use is recommended but not required. The **goodHearing A1** comes pre-charged however initial charge will vary.



- 2) 3+ hours for full charge gives up to 10 hours of usage time



3) When plugged into a power source



Solid RED LED Indicator : Charging



Solid BLUE LED Indicator : Fully charged

3.3) Switching Between Modes

1) There are 3 pre-set modes on the **goodHearing A1**.

- General : Standard Mode


Mostly used indoors where the environment is not too noisy.

- Conversation : Noise Cancelling Mode

Used in louder environment such as outdoor public areas, restaurants.

- In-car : High-level Noise Cancelling Mode

Providing superior noise cancelling in even the loudest environments, such as in transportation.

2) To switch between the 3 modes click the  button (located on the center of the **goodHearing A1** device). A beep will be heard when switching between the 3 modes to indicate which mode it is currently set.

- Mode 1 : A single beep will be heard
(General: Standard Mode)
- Mode 2 : Two beeps will be heard
(Conversation: Noise Cancelling Mode)
- Mode 3 : Three beeps will be heard
(In-car: High-level Noise Cancelling Mode)

4. PAIRING BLUETOOTH

4.1) Download App

Please scan QR code to download the App or search “**goodHearing**” on App Store / Google Play



goodHearing




iOS




Android

4.2) Bluetooth Connect


1. Turn on **goodHearing A1** and have both devices within arm length of each other.
2. Press and hold the **(M)** for 3 seconds, the LED Indicator will flash alternating BLUE and RED colors.
3. Turn on the Bluetooth function on your mobile device.



4. Scan for nearby devices, search and pair with **goodHearing A1**.



5. When pairing has been successful your phone will show that the status is “Connected” with **goodHearing A1**.



- ✘ The pairing steps may vary according to your mobile operating system.

4.3) Notes

- Paired devices will remember each other unless pairing information is deleted manually, so there is no need to pair again. The **goodHearing A1** can pair and remember 2 Bluetooth devices.
- Repeat the steps outlined above to set-up your **goodHearing A1** with another Bluetooth device.
- When a 3rd Bluetooth device is paired, the oldest pairing device information with **goodHearing A1** will be deleted.
- The **goodHearing A1** will connect automatically with the latest connected Bluetooth device.
- If you want to use another Bluetooth device to pair with **goodHearing A1**, switch off the Bluetooth function of the current Bluetooth device, connect the new Bluetooth mobile device to the unit by selecting **goodHearing A1** on the Bluetooth device list and press connect.

- When the Bluetooth device that is connected to the **goodHearing A1** is not playing any sort of media, the device will automatically switch back to being its primary function until the connected device starts playing media or receives a phone call again.
- Should the Bluetooth be disconnected due to range or any other reason, reconnect the **goodHearing A1** with the device that you wish to pair it with.
- If you have connected the **goodHearing A1** to a device and the 2 devices gets out of connecting range, the **goodHearing A1** will have a reminder tone in the background. To switch this tone off, simple turn off the unit for 5 seconds and turn it back on.


5. SAFETY & MAINTENANCE

5.1) Earbud Maintenance

- Clean your ear thoroughly before each use.



- Remove all remaining debris from the earbud and ensure that the hole on the side is also free from any debris after each use.




- ※ Use the Cleaning Brush to clean any earwax adhered on the hearing aid, especially in and around the hearing aid and earbud tip.








5.2) Battery Maintenance

- Do not consistently discharge the battery completely before recharging as this may cause depletion of the batteries life.
- We suggest recharging **goodHearing A1** regardless of its battery status before sleeping at night.
- Do not attempt to remove the battery as this may cause damage to the **goodHearing A1** device or may cause injuries.
- Ensure the device is fully charged before removing it from charge.

5.3) LED Indicator status

- Solid RED: Charging
- Solid BLUE: Fully Charged
- Flashing BLUE: ON (working normally)
- Flashing RED and BLUE:
Searching for Bluetooth devices to be paired
- Flashing PURPLE:
Bluetooth Connected with media playing



LED Indicator	
	Red :  ON / OFF
	Blue : Bluetooth Standby
	Red & Blue : Bluetooth Search
	Red-Purple-Blue : Bluetooth Connected with media playing
	Blue : Fully Charged
	Red : Charging

6. FAQ

No Sound ?

- Turn on the device
- Increase the volume
- Clean out the earphone caps
- Check if the earphone cable damages
- Recharge the device

Charging problems?

- Check that the USB is connected correctly
- Ensure that the power supply is normal
(wall outlet switch is on)

Is it compatible for both Android and iOS systems ?

Yes, is compatible for both Android and iOS systems

Can't power on

Out of battery and needs to be charged

There is a lag when I change from listening to media on my phone and back to the hearing aid function is this normal?

Yes, there may be up to 5 second delay depends on the operating system you currently running on your phone

Can I pick up phone calls with the hearing aid ?

Yes, you can, depend on Android or iOS system, your phones will have a default tone that the hearing aid will play as a ringtone

There is still no sound on my SONY phone ?

Please turn off the following sound settings:

Screen lock sound, Vibration on touch,

Keyboard touch tone

There is still no sound on my HTC phone ?

Please turn off the following sound settings:

Touch sounds, Screen lock sounds, Pull-to-

refresh sound, Keyboard touch tones

7. SPECIFICATION

Product name goodHearing A1

Model name UP-6EB3

goodHearing A1	
Dimensions	160 x 190 x 11 mm
Weight	35 g / pcs
Color	White / Black
Battery Type	Rechargeable Battery
Continuous Use	Up to 10 hours per charge
Listening Mode	3 Modes
Volume	10 Levels
Storage Temperature	14oF ~ 140oF
Storage Humidity	30% ~ 85%

8. WARNING AND REMINDER

- Keep the device dry
- Keep away from direct sunlight
- Keep away from reach of children
- Keep away from high temperature areas
- Avoid liquid near the device
- Avoid discharging the battery for long periods
- Avoid overcharging the device for excessive time periods
- Do not use the device at high volumes
- Do not listen for excessive time periods
- Do not attempt to fix the unit yourself, take it to an authorized dealer
- Ensure the earphones are fitted securely in the ear to reduce sound leak

9. LEGAL INFORMATION

Warranty

1. Mimitakara provides a one-year warranty,
2. Under normal conditions use, provides free repair services. Earplugs, horns, horn lines, (Personal hygiene accessories are not covered by the warranty)
3. The product improperly handled/ miss-used/ damaged shall be charged at the discretion when they need to be repaired.
4. Serious damage product or self-disassembly or modification will not be repaired.
5. Damage caused by uncontrolled factors such as natural disasters and fires is not covered in the warranty.
6. If you need repair services, please send the product with the warranty form to the company.
7. 30 day risk free money back guarantee
8. 31 to 60 days restocking fee

FCC

Federal Communications Commission (FCC) Statement

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

FCC RF Radiation Exposure Statement:


1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

2. For portable operation, this device has been tested and meets FCC RF exposure guidelines. When used with an accessory that contains metal may not ensure compliance with FCC RF exposure guidelines.

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- 1) this device may not cause harmful interference and
- 2) this device must accept any interference received, including interference that may cause undesired operation.



WARNING: If you are younger than 18, do not use this.

You should go to a doctor, preferably an ear-nose-throat doctor (an ENT), because your condition needs specialized care.

Over-the-counter hearing aids are only for users who are age 18 or older.

This hearing aid is for adults with signs of mild to moderate hearing loss. How do you know if you have this?

- You have trouble hearing speech in noisy places
- You find it hard to follow speech in groups
- You have trouble hearing on the phone
- Listening makes you tired
- You need to turn up the volume on the TV or radio, and other people complain it's too loud

Some people with hearing loss may need help from a hearing healthcare professional. How do you know if you need to see one?

- You can't hear speech even if the room is quiet
- You don't hear loud sounds well, for example, you don't hear loud music, power tools, engines, or other very noisy things.

If your hearing loss makes it hard to hear loud noises, this hearing aid may not be your best choice without help from a professional. If this hearing aid does not help you enough, ask for help from a hearing healthcare professional.

WARNING: When to See a Doctor

If you have any of the problems listed below, please see a doctor, preferably an ear-nose-throat doctor (an ENT).

- Your ear has a birth defect or an unusual shape.

- Your ear was injured or deformed in an accident.
- You saw blood, pus, or fluid coming out of your ear in the past 6 months
- Your ear feels painful or uncomfortable
- You have a lot of ear wax, or you think something could be in your ear
- You get really dizzy or have a feeling of spinning or swaying (called vertigo)
- Your hearing changed suddenly in the past 6 months
- Your hearing changes: it gets worse then gets better again
- You have worse hearing in one ear
- You hear ringing or buzzing in only one ear

This information and other labeling, including the user instructional brochure, are available on the internet at: www.mimitakarahealth.com

You may also call +1 949 695 7681 (9am-6pm PST) or write to info@mimitakara.com or to request a paper copy of this information and other labeling.

Manufacturer's return policy:

- 30-Day Money-Back Guarantee
- 1-Year Manufacturer's Warranty

See return details on the internet at www.mimitakarahealth.com



FCC ID: 2BBAZ-6EBX

The device complies with Part 15 of the FCC. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.